**EN** EXPERT CLEANING, SPECIALLY DESIGNED FOR ORTHODONTIC APPLIANCES WEARERS. PROVEN SUPERIOR EFFICACY:

Precision bristles provide 45% more efficacy compared to a manual orthodontic toothbrush<sup>1</sup>: ultra-thin tipped, longer bristles effectively clean in hard-to-reach spaces, end-rounded shorter bristles allow an effective clean on teeth surfaces and around brackets. SONIC vibrations (12500/min) help remove plaque. Portable: light, cable-free, cap included. Long lasting battery – Up to 90 days of use with 1 AAA battery<sup>2</sup>.

Scan our QR code or visit SunstarGUM.com to learn how to take the best care of your orthodontic appliances.

**(A) HOW TO REPLACE THE BRUSH HEAD:** Holding base, pull on head. Place new head on base and push down. Compatible with all GUM<sup>®</sup> SONIC refill heads.

**(B) HOW TO REPLACE THE BATTERY:** Turn and pull-down handle base. Replace it with one new AAA alkaline battery only. When replacing the battery, align symbols (a and b) on the bottom part for correct positioning. Close the container and dispose the battery according to local regulations.

**HOW TO USE THE TOOTHBRUSH:** Apply toothpaste on the brush head. Insert the brush head inside the mouth. Turn on the brush by pressing the upper central button. Use the same brushing technique as with a manual toothbrush without applying any additional force. Make sure you also carefully brush over and below the braces/bands. Press the lower central button to turn it on. Ask additional advice to your orthodontist or dentist on how best to brush with braces.

It is recommended to change the brush head every 3 months, when bristles are worn out or after sickness.

**HOW TO CLEAN THE TOOTHBRUSH:** Turn off the brush. Rinse brush head with water after brushing.

**WARNING.** Keep out of reach of children. Choking hazard. Small parts. Not for children under 3 years. For children 3 years and older, use it under adult supervision only.

**PACKAGING CONTENT:** 1 toothbrush handle, 1 head, 1 AAA battery, 1 toothbrush cap.

<sup>1</sup> Tested *in vitro*  
<sup>2</sup> with 2 min brushing session twice per day.

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Jakelija: Sunstar Sverige AB, SE-436 33 Askim
Dystrybuowana przez: Sunstar Europe S.A. Oddział w Polsce, PL- 02-222, Warszawa
Distribuerad av: Sunstar Sverige AB, SE-436 33 Askim
Vertrieb durch: Sunstar Deutschland GmbH, DE-79677 Schönau